Membracel

Porous Regenerative Membrane

Membrana Regeneradora Porosa

Pediatric and adult use

Sterile — sterilized using ethylene oxide

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Description • Descripción

Membracel – Porous Regenerative Membrane consists of interwoven crystalline cellulose microfibrils, it has artificially developed pores, fine texture and uniform structure. These features allow protection from skin lesions and promote tissue regeneration by the acceleration of granulation tissue formation and subsequent epithelialization.

Product expected action: relief from pain, injury protection, drainage of secretions, to reduce the exposed region to infectious agents and accelerating the healing process.

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Indications • Indicaciones

The **Porous Regenerative Membrane** works as a temporary replacement of the skin and is suitable for injury treatments resulting from the loss of epithelium that are characterized as superficial or deep wound, with little or abundant exudation, such as second-degree burns, dermal abrasions, excoriation, recipient and donor areas for skin grafts, nail beds (post-nail excision), pressure sores, varicose ulcers of stasis, arterial ulcers, plantar perforating maladies, epidermolysis bullosa, post-physical cauterization (cryotherapy, thermal cauterization, CO₂ laser and Erbium).

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Contraindications • Contraindicaciones

The **Porous Regenerative Membrane** should not be used to control bleeding, tunnel lesions, malignant lesions or suspected malignancy.

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Choosing the membrane • Escogiendo la membrana

Choosing the right **Porous Regenerative Membrane** is essential for treatment.

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Consider the size and shape of the wound to choose the most **suitable shape and size**, so that the membrane covers the entire exposed injury area and exceeding by at least 1 cm its edges:

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Consider the type of wound, the volume and characteristics of the existing exudate in the wound, to choose the right pore size:

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AVERAGE – 1 to 2 mm
Suitable for low and/or moderate exudation. Used in most cases.

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LARGE – 2 to 3 mm
Indicated in cases of abundant and/or thick exudation.

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The artificially produced pores in **Membracel – Porous Regenerative Membrane** have the function of allowing the free passage of exudate from the lesion bed to the secondary dressing, which must have absorbent characteristics.

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Method of use • Modo de usar

1 Proceed to clean the lesion with saline solution (0.9%) in jet or another specific solution for cleaning; if necessary, debridement to remove non-viable tissue.

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2 Remove the **Porous Regenerative Membrane** from envelope and place it on the wound so as to cover the lesion with an excess of 1 cm beyond its edges. The product can be applied on both sides.

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3 With the use of forceps and sterile gauze, perfectly settle the product over the injured area, gently pressing the gauze on the membrane.

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4 Wait 3 to 4 minutes for the occurrence of membrane self-adherence in the lesion and the peripheral region.

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5 If necessary, apply an absorbent secondary dressing on the membrane, such as sterile gauze. For fixation, use bands, compression bandages and microporous tape.

For greater efficiency in the treatment, it is recommended to change the secondary dressing whenever it is saturated with exudate or lose the ability to absorb, or in accordance with the guidance of health professionals.

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Intervals for membrane change • Intervalos de cambio de la membrana

1 The change frequency will depend on the condition of the wound and the volume of exudate. Initially it is normal that the need for change will occur every 5 days. Over time, this interval will naturally increase, due to the reduction in the exudate volume secreted by the injury.

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2 The porous membrane can be maintained for up to 12 days if it does not occur peel off, disruption, signs of infection, acute pain, foul odor or stench.

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Procedures for changing secondary dressing •

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- 1 To remove secondary dressing applied over the membrane.

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- 2 Cleaning the wound area without removing the membrane using saline solution (0.9%) and a gauze to compress and remove the excess exudate under the membrane, through its pores. In this procedure, avoid too moisten the membrane to prevent its detachment.

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- 3 Place gauze or other sterile absorbent secondary dressing. Secure the secondary dressing with microporous adhesive tape or band, so that it is lightly pressed over the membrane/lesion.

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Procedures for changing porous membrane •

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If you need to change the **Porous Regenerative Membrane**, adopt the following:

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- 1 Moisten the adhered part of the membrane with gauze soaked in saline solution (0.9%) for about 4 to 5 minutes.

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- 2 Proceed to peel off the edge of the membrane, with gentle rubbing of gauze soaked in saline solution (0.9%) in the periphery-center direction of the injury until complete removal of the membrane.

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- 3 Clean the wound bed using saline solution (0.9%) in jet or another specific solution. Apply a new membrane, observing the procedures described in Method of use.

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Special care • Cuidados especiales

- 1 Never use oily products, ointments, creams or gels under the membrane, since this practice prevents its attachment to the lesion. In need for additional topical treatment, the products in the form of aqueous solution are the most recommendable. However, avoid the use of hydrogen peroxide, as it will cause chemical degradation of the membrane.

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- 2 Avoid wetting the membrane during bath (protect the dressing with waterproof material).

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- 3 **Attention:** sterilization is guaranteed as long as the packaging is not damaged or opened prior to use. Do not use the **Porous Regenerative Membrane** if the package is broken off and/or out-dated.

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- 4 Systemic antimicrobial therapy can be prescribed by your doctor when the wound presents infection. In this case, the **Porous Regenerative Membrane** can be used in conjunction with systemic antibiotics.

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- 5 If necessary to cut the **Porous Regenerative Membrane** dressing, it's recommended to use sterile scissors.

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- 6 The use of **Porous Regenerative Membrane** should be supervised by a health professional.

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Storage • Almacenamiento

This product, when stored at room temperature (between 15°C/ 59°F and 30°C/ 86°F), away from light and moisture, has a validity of 60 months from the date of manufacture.

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Packs with 1 or 10 **Porous Regenerative Membranes**.

Envase con 1 o 10 **Membranas Regeneradoras Porosas**.

SYMBOL	REFERENCE
	Expiration date Fecha de vencimiento
	Manufactured by Fabricante
	Lot number Código de lote
	Sterile product Producto estéril
	Sterilized using ethylene oxide Esterilizado con óxido de etileno
	Product for single use only Producto de uso único
	Attention, see instructions for use Precaución, consultar los documentos adjuntos
	Store at temperature between 59°F (15°C) and 86°F (30°C) Límite de temperatura entre 59°F (15°C) y 86°F (30°C)
	Store in a dry place Almacene en lugar seco
	Do not use if package is damaged No utilizar si el envase está dañado

STERILIZED USING ETHYLENE OXIDE. PRODUCT FOR SINGLE USE ONLY - DO NOT REUSE. DO NOT USE IF PACKAGE IS OPEN OR DAMAGE. FOLLOW THE DIRECTIONS OF YOUR PHYSICIAN.

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Lot number, manufacturing date and expiration date:
see box or envelope.

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VUELO

Manufactured and distributed by:
Fabricado y distribuido por:

Vuelo Pharma Ltda.

Rua Antônio de Paula Braga, 83
Almirante Tamandaré - Paraná
CEP: 83506-010
CNPJ: 04.159.536/0001-05

vuelopharma.com

sac@vuelopharma.com
Phone/Tel.: +55 41 3657-7611
SAC: 41 3657-9342

Technical support

Responsable técnico
Andréa P. Canuto von der Heyde
CRF-PR 18.418

ANVISA Registration
80074780002

Approved
3006611176